

# Module Interface Specification for Software Engineering

Team 17, Team RAdiAIdance

Allison Cook

Ibrahim Issa

Mohaansh Pranjali

Nathaniel Hu

Tushar Aggarwal

January 17, 2024

# 1 Revision History

Date	Version	Notes
01/11/2024	0.0	Initial Document
01/12/2024	0.1	Started adding MIS for Module sections
01/14/2024	0.2	Added in Table from Module Guide (MG) into Module Decomposition section
01/16/2024	0.3	Made small updates to Table 1 in Module Decomposition section; Continued adding MIS for Module sections
01/17/2024	0.4	Completed the Symbols, Abbreviations and Acronyms and Introduction sections; Finished adding in MIS for Module sections

## 2 Symbols, Abbreviations and Acronyms

This section records the symbols, abbreviations and acronyms information for easy reference for terms used in this document.

For information on most of the symbols, abbreviations and acronyms referenced in this document, see the SRS Documentation at the following link: <https://github.com/tusharagg1/chest-x-ray-ai/blob/main/docs/SRS/SRS.pdf>.

The information on the rest of the symbols, abbreviations and acronyms referenced in this document are shown in the table below.

symbol	description
AI/ML	Artificial Intelligence/Machine Learning
DICOM	Digital Imaging and Communications in Medicine; technical standard for digital storage/transmission of medical images and related information
GUI	Graphical User Interface
JPEG/JPG	Joint Photographic Experts Group; digital image compression standard, image format
M	Module
MG	Module Guide
MVC	Model-View-Controller Software Architecture
NLP	Natural Language Processing
SRS	Software Requirements Specification
Software Engineering	The Process of Designing and Developing Software; a reference to the software application described in this document

# Contents

<b>1</b>	<b>Revision History</b>	<b>i</b>
<b>2</b>	<b>Symbols, Abbreviations and Acronyms</b>	<b>ii</b>
<b>3</b>	<b>Introduction</b>	<b>1</b>
<b>4</b>	<b>Notation</b>	<b>1</b>
<b>5</b>	<b>Module Decomposition</b>	<b>1</b>
<b>6</b>	<b>MIS of Medical Institution Interface Module</b>	<b>3</b>
6.1	Module . . . . .	3
6.2	Uses . . . . .	3
6.3	Syntax . . . . .	3
6.3.1	Exported Constants . . . . .	3
6.3.2	Exported Access Programs . . . . .	3
6.4	Semantics . . . . .	3
6.4.1	State Variables . . . . .	3
6.4.2	Environment Variables . . . . .	3
6.4.3	Assumptions . . . . .	3
6.4.4	Access Routine Semantics . . . . .	3
6.4.5	Local Functions . . . . .	4
<b>7</b>	<b>MIS of Chest X-Ray Read Module</b>	<b>5</b>
7.1	Module . . . . .	5
7.2	Uses . . . . .	5
7.3	Syntax . . . . .	5
7.3.1	Exported Constants . . . . .	5
7.3.2	Exported Access Programs . . . . .	5
7.4	Semantics . . . . .	5
7.4.1	State Variables . . . . .	5
7.4.2	Environment Variables . . . . .	5
7.4.3	Assumptions . . . . .	5
7.4.4	Access Routine Semantics . . . . .	5
7.4.5	Local Functions . . . . .	6
<b>8</b>	<b>MIS of Results Generation Module</b>	<b>7</b>
8.1	Module . . . . .	7
8.2	Uses . . . . .	7
8.3	Syntax . . . . .	7
8.3.1	Exported Constants . . . . .	7
8.3.2	Exported Access Programs . . . . .	7

8.4	Semantics . . . . .	7
8.4.1	State Variables . . . . .	7
8.4.2	Environment Variables . . . . .	7
8.4.3	Assumptions . . . . .	7
8.4.4	Access Routine Semantics . . . . .	7
8.4.5	Local Functions . . . . .	8
<b>9</b>	<b>MIS of Report Component Generation Module</b>	<b>9</b>
9.1	Module . . . . .	9
9.2	Uses . . . . .	9
9.3	Syntax . . . . .	9
9.3.1	Exported Constants . . . . .	9
9.3.2	Exported Access Programs . . . . .	9
9.4	Semantics . . . . .	9
9.4.1	State Variables . . . . .	9
9.4.2	Environment Variables . . . . .	9
9.4.3	Assumptions . . . . .	9
9.4.4	Access Routine Semantics . . . . .	9
9.4.5	Local Functions . . . . .	10
<b>10</b>	<b>MIS of Database Operations Module</b>	<b>11</b>
10.1	Module . . . . .	11
10.2	Uses . . . . .	11
10.3	Syntax . . . . .	11
10.3.1	Exported Constants . . . . .	11
10.3.2	Exported Access Programs . . . . .	11
10.4	Semantics . . . . .	11
10.4.1	State Variables . . . . .	11
10.4.2	Environment Variables . . . . .	11
10.4.3	Assumptions . . . . .	11
10.4.4	Access Routine Semantics . . . . .	11
10.4.5	Local Functions . . . . .	12
<b>11</b>	<b>MIS of User Authentication/Management Module</b>	<b>13</b>
11.1	Module . . . . .	13
11.2	Uses . . . . .	13
11.3	Syntax . . . . .	13
11.3.1	Exported Constants . . . . .	13
11.3.2	Exported Access Programs . . . . .	13
11.4	Semantics . . . . .	13
11.4.1	State Variables . . . . .	13
11.4.2	Environment Variables . . . . .	13
11.4.3	Assumptions . . . . .	13

11.4.4	Access Routine Semantics . . . . .	14
11.4.5	Local Functions . . . . .	15
<b>12</b>	<b>MIS of App GUI Module</b>	<b>16</b>
12.1	Module . . . . .	16
12.2	Uses . . . . .	16
12.3	Syntax . . . . .	16
12.3.1	Exported Constants . . . . .	16
12.3.2	Exported Access Programs . . . . .	16
12.4	Semantics . . . . .	16
12.4.1	State Variables . . . . .	16
12.4.2	Environment Variables . . . . .	16
12.4.3	Assumptions . . . . .	16
12.4.4	Access Routine Semantics . . . . .	17
12.4.5	Local Functions . . . . .	17
<b>13</b>	<b>MIS of Login Module</b>	<b>18</b>
13.1	Module . . . . .	18
13.2	Uses . . . . .	18
13.3	Syntax . . . . .	18
13.3.1	Exported Constants . . . . .	18
13.3.2	Exported Access Programs . . . . .	18
13.4	Semantics . . . . .	18
13.4.1	State Variables . . . . .	18
13.4.2	Environment Variables . . . . .	18
13.4.3	Assumptions . . . . .	18
13.4.4	Access Routine Semantics . . . . .	18
13.4.5	Local Functions . . . . .	19
<b>14</b>	<b>MIS of Perform Scan Module</b>	<b>20</b>
14.1	Module . . . . .	20
14.2	Uses . . . . .	20
14.3	Syntax . . . . .	20
14.3.1	Exported Constants . . . . .	20
14.3.2	Exported Access Programs . . . . .	20
14.4	Semantics . . . . .	20
14.4.1	State Variables . . . . .	20
14.4.2	Environment Variables . . . . .	20
14.4.3	Assumptions . . . . .	20
14.4.4	Access Routine Semantics . . . . .	20
14.4.5	Local Functions . . . . .	21

<b>15 MIS of View Results Module</b>	<b>22</b>
15.1 Module . . . . .	22
15.2 Uses . . . . .	22
15.3 Syntax . . . . .	22
15.3.1 Exported Constants . . . . .	22
15.3.2 Exported Access Programs . . . . .	22
15.4 Semantics . . . . .	22
15.4.1 State Variables . . . . .	22
15.4.2 Environment Variables . . . . .	22
15.4.3 Assumptions . . . . .	22
15.4.4 Access Routine Semantics . . . . .	22
15.4.5 Local Functions . . . . .	23
<b>16 MIS of AI Model Module</b>	<b>24</b>
16.1 Module . . . . .	24
16.2 Uses . . . . .	24
16.3 Syntax . . . . .	24
16.3.1 Exported Constants . . . . .	24
16.3.2 Exported Access Programs . . . . .	24
16.4 Semantics . . . . .	24
16.4.1 State Variables . . . . .	24
16.4.2 Environment Variables . . . . .	24
16.4.3 Assumptions . . . . .	24
16.4.4 Access Routine Semantics . . . . .	25
16.4.5 Local Functions . . . . .	25
<b>17 MIS of NLP Model Module</b>	<b>26</b>
17.1 Module . . . . .	26
17.2 Uses . . . . .	26
17.3 Syntax . . . . .	26
17.3.1 Exported Constants . . . . .	26
17.3.2 Exported Access Programs . . . . .	26
17.4 Semantics . . . . .	26
17.4.1 State Variables . . . . .	26
17.4.2 Environment Variables . . . . .	26
17.4.3 Assumptions . . . . .	26
17.4.4 Access Routine Semantics . . . . .	26
17.4.5 Local Functions . . . . .	27
<b>18 MIS of Backend Module</b>	<b>28</b>
18.1 Module . . . . .	28
18.2 Uses . . . . .	28
18.3 Syntax . . . . .	28

18.3.1	Exported Constants . . . . .	28
18.3.2	Exported Access Programs . . . . .	28
18.4	Semantics . . . . .	28
18.4.1	State Variables . . . . .	28
18.4.2	Environment Variables . . . . .	28
18.4.3	Assumptions . . . . .	28
18.4.4	Access Routine Semantics . . . . .	28
18.4.5	Local Functions . . . . .	29
<b>19</b>	<b>MIS of App Controller Module</b>	<b>30</b>
19.1	Module . . . . .	30
19.2	Uses . . . . .	30
19.3	Syntax . . . . .	30
19.3.1	Exported Constants . . . . .	30
19.3.2	Exported Access Programs . . . . .	30
19.4	Semantics . . . . .	30
19.4.1	State Variables . . . . .	30
19.4.2	Environment Variables . . . . .	30
19.4.3	Assumptions . . . . .	30
19.4.4	Access Routine Semantics . . . . .	31
19.4.5	Local Functions . . . . .	31
<b>20</b>	<b>Appendix</b>	<b>33</b>



### 3 Introduction

The following document details the Module Interface Specifications for the [Your Program Name Here] software application. This software application (sometimes referred to as Software Engineering in this document) performs scans of chest x-ray images, looking for diseases/infections and making predictions. Those scan results and predictions of diseases/infections are then put into natural language radiology reports (or components) and returned.

Complementary documents include the System Requirement Specifications and Module Guide. The full documentation and implementation can be found at <https://github.com/tusharagg1/chest-x-ray-ai/tree/main>.

### 4 Notation

The structure of the MIS for modules comes from Hoffman and Strooper (1995), with the addition that template modules have been adapted from Ghezzi et al. (2003). The mathematical notation comes from Chapter 3 of Hoffman and Strooper (1995). For instance, the symbol  $:=$  is used for a multiple assignment statement and conditional rules follow the form  $(c_1 \Rightarrow r_1 | c_2 \Rightarrow r_2 | \dots | c_n \Rightarrow r_n)$ .

The following table summarizes the primitive data types used by Software Engineering.

Data Type	Notation	Description
character	char	a single symbol or digit
integer	$\mathbb{Z}$	a number without a fractional component in $(-\infty, \infty)$
natural number	$\mathbb{N}$	a number without a fractional component in $[1, \infty)$
real	$\mathbb{R}$	any number in $(-\infty, \infty)$

The specification of Software Engineering uses some derived data types: sequences, strings, and tuples. Sequences are lists filled with elements of the same data type. Strings are sequences of characters. Tuples contain a list of values, potentially of different types. In addition, Software Engineering uses functions, which are defined by the data types of their inputs and outputs. Local functions are described by giving their type signature followed by their specification.

### 5 Module Decomposition

The following table is taken directly from the Module Guide document for this project.

Level 1	Level 2
Hardware-Hiding Module	MedInstInter
Behaviour-Hiding Module	ChestXRayRead
	ResultsGen
	RepCompGen
	DatabaseOps
	UserAuthMgmt
	Login
	PerfScan
Software Decision Module	ViewResults
	AIModel
	NLPModel
	Backend
	AppController
	AppGUI

Table 1: Module Hierarchy

## 6 MIS of Medical Institution Interface Module

### 6.1 Module

MedInstInter

### 6.2 Uses

N/A

### 6.3 Syntax

#### 6.3.1 Exported Constants

N/A

#### 6.3.2 Exported Access Programs

Name	In	Out	Exceptions
connectToInst	instID: str, creden- tials: str	connectionStatus: bool	InvalidCredentialsException, InstNot- FoundException

### 6.4 Semantics

#### 6.4.1 State Variables

- connectedInsts: Set(str) - maintains a set of connected institution IDs.

#### 6.4.2 Environment Variables

- InstsITSys: Set(str) - the set of external IT systems the application interfaces with to retrieve/exchange information.

#### 6.4.3 Assumptions

- Patient data is stored in the medical institution's database, and the software intends to interface with their server to access that information.

#### 6.4.4 Access Routine Semantics

connectToInst():

- transition:
  - Adds 'instID' to 'connectedInsts' if the provided 'credentials' is valid.

- output:
  - ‘connectionStatus’ is set to True if the connection is successful, False otherwise.
- exception:
  - Throws ‘InvalidCredentialsException’ if the provided credentials are invalid.
  - Throws ‘InstNotFoundException’ if the specified ‘instID’ does not exist.

#### **6.4.5 Local Functions**

N/A

## 7 MIS of Chest X-Ray Read Module

### 7.1 Module

ChestXRayRead

### 7.2 Uses

N/A

### 7.3 Syntax

#### 7.3.1 Exported Constants

N/A

#### 7.3.2 Exported Access Programs

Name	In	Out	Exceptions
processImg	img: JPEG image	procImg: Image	Processed InvalidImageFormatException

### 7.4 Semantics

#### 7.4.1 State Variables

N/A

#### 7.4.2 Environment Variables

N/A

#### 7.4.3 Assumptions

N/A

#### 7.4.4 Access Routine Semantics

processImg():

- transition:
  - Initiates the AIModel module to process the provided ‘img’.
- output:

- ‘procImg’ contains the processed information and findings from the chest X-ray analysis.
- exception:
  - Throws ‘InvalidImageFormatException’ if the provided ‘image’ is not a JPEG or JPG image.

#### **7.4.5 Local Functions**

N/A

## 8 MIS of Results Generation Module

### 8.1 Module

ResultsGen

### 8.2 Uses

N/A

### 8.3 Syntax

#### 8.3.1 Exported Constants

N/A

#### 8.3.2 Exported Access Programs

Name	In	Out	Exceptions
generateResults	procImg: Image	Processed classification: Disease Classification	-

### 8.4 Semantics

#### 8.4.1 State Variables

N/A

#### 8.4.2 Environment Variables

N/A

#### 8.4.3 Assumptions

N/A

#### 8.4.4 Access Routine Semantics

generateResults():

- transition:
  - Utilizes the AIModel module to interpret the processed image and generate a disease classification.
- output:

- ‘classification’ contains the generated disease classification for each disease.
- exception: N/A

#### **8.4.5 Local Functions**

N/A



## 9 MIS of Report Component Generation Module

### 9.1 Module

RepCompGen

### 9.2 Uses

N/A

### 9.3 Syntax

#### 9.3.1 Exported Constants

N/A

#### 9.3.2 Exported Access Programs

Name	In	Out	Exceptions
generateReport	diagnosis: Disease Diagnosis	report: Radiology Report	-

### 9.4 Semantics

#### 9.4.1 State Variables

N/A

#### 9.4.2 Environment Variables

N/A

#### 9.4.3 Assumptions

#### 9.4.4 Access Routine Semantics

generateReport():

- transition: N/A
- output:
  - ‘report’ contains the generated radiology report based on the provided disease diagnosis.
- exception: N/A

#### 9.4.5 Local Functions

N/A

## 10 MIS of Database Operations Module

### 10.1 Module

DatabaseOps

### 10.2 Uses

N/A

### 10.3 Syntax

#### 10.3.1 Exported Constants

N/A

#### 10.3.2 Exported Access Programs

Name	In	Out	Exceptions
storeReport	report: Radiology Report, patientID: str	success: bool	ReportStorageException
retrieveReport	patientID: str	report: Radiology Re- port	ReportRetrievalException

### 10.4 Semantics

#### 10.4.1 State Variables

- ‘connectDatabase: bool’ indicates whether the module is currently connected to the database.

#### 10.4.2 Environment Variables

N/A

#### 10.4.3 Assumptions

N/A

#### 10.4.4 Access Routine Semantics

storeReport():

- transition:

- Stores the provided ‘report’ in the database associated with the specified ‘patientID’.

- output:

- ‘success’ is set to True if the storing operation is successful, False otherwise.

- exception:

- Throws ‘ReportStorageException’ if there is an issue storing the report.

retrieveReport():

- transition:

- Retrieves the radiology report associated with the specified ‘patientID’ from the database.

- output:

- ‘report’ contains the retrieved radiology report.

- exception:

- Throws ‘ReportRetrievalException’ if there is an issue retrieving the report.

#### 10.4.5 Local Functions

N/A

## 11 MIS of User Authentication/Management Module

### 11.1 Module

UserAuthMgmt

### 11.2 Uses

N/A

### 11.3 Syntax

#### 11.3.1 Exported Constants

N/A

#### 11.3.2 Exported Access Programs

Name	In	Out	Exceptions
authenticateUser	username: str, password: str	status: bool	InvalidCredentialsException, UserNotFoundException
createUserAccount	username: str, password: str	success: bool	UserCreationException
deleteUserAccount	username: str, password: str	success: bool	UserDeletionException
checkAuthentication	username: str	isAuthorized: bool	-

### 11.4 Semantics

#### 11.4.1 State Variables

N/A

#### 11.4.2 Environment Variables

N/A

#### 11.4.3 Assumptions

N/A

#### 11.4.4 Access Routine Semantics

authenticateUser():

- transition:
  - Verifies the provided ‘username’ and ‘password’ for authentication.
- output:
  - ‘status’ is set to True if authentication is successful, False otherwise.
- exception:
  - Throws ‘InvalidCredentialsException’ if the provided credentials are invalid.
  - Throws ‘UserNotFoundException’ if the specified user is not found.

createUserAccount():

- transition:
  - Creates a user account with the provided ‘username’ and ‘password’.
- output:
  - ‘success’ is set to True if the account creation is successful, False otherwise.
- exception:
  - Throws ‘UserCreationException’ if there is an issue creating the user account.

deleteUserAccount():

- transition:
  - Deletes the user account associated with the specified ‘username’.
- output:
  - ‘success’ is set to True if the account deletion is successful, False otherwise.
- exception:
  - Throws ‘UserDeletionException’ if there is an issue deleting the user account.

checkAuthentication():

- transition:
  - Checks whether the specified ‘username’ is currently authorized.

- output:
  - ‘isAuthorized’ is set to True if the user is authorized, False otherwise.
- exception: N/A

#### **11.4.5 Local Functions**

N/A

## 12 MIS of App GUI Module

### 12.1 Module

AppGUI

### 12.2 Uses

- Login
- PerfScan
- ViewResults

### 12.3 Syntax

#### 12.3.1 Exported Constants

N/A

#### 12.3.2 Exported Access Programs

Name	In	Out	Exceptions
displayLoginPage	-	-	-
displayScanPage	-	-	-
displayResultsPage	-	-	-

### 12.4 Semantics

#### 12.4.1 State Variables

N/A

#### 12.4.2 Environment Variables

N/A

#### 12.4.3 Assumptions

N/A



#### 12.4.4 Access Routine Semantics

displayLoginPage():

- transition: Navigates to and displays the login page for the application.
- output: N/A
- exception: N/A

displayScanPage():

- transition: Navigates to and displays the page for inputting an x-ray image for scanning.
- output: N/A
- exception: N/A

displayResultsPage():

- transition: Navigates to and displays the page for viewing scan results and reports.
- output: N/A
- exception: N/A

#### 12.4.5 Local Functions

N/A

## 13 MIS of Login Module

### 13.1 Module

Login

### 13.2 Uses

N/A

### 13.3 Syntax

#### 13.3.1 Exported Constants

N/A

#### 13.3.2 Exported Access Programs

Name	In	Out	Exceptions
login	username: str, password: str	loginStatus: bool	InvalidCredentialsException, UserNotFoundException

### 13.4 Semantics

#### 13.4.1 State Variables

N/A

#### 13.4.2 Environment Variables

N/A

#### 13.4.3 Assumptions

N/A

#### 13.4.4 Access Routine Semantics

login():

- transition:
  - Authenticates the provided 'username' and 'password'.
- output:

- 'loginStatus' is set to True if login is successful, False otherwise.
- exception:
  - Throws 'InvalidCredentialsException' if the provided credentials are invalid.
  - Throws 'UserNotFoundException' if the specified user is not found.

#### **13.4.5 Local Functions**

N/A

## 14 MIS of Perform Scan Module

### 14.1 Module

PerfScan

### 14.2 Uses

N/A

### 14.3 Syntax

#### 14.3.1 Exported Constants

N/A

#### 14.3.2 Exported Access Programs

Name	In	Out	Exceptions
initiateScan	img: X-Ray Image	-	InvalidImageFormatException

### 14.4 Semantics

#### 14.4.1 State Variables

N/A

#### 14.4.2 Environment Variables

N/A

#### 14.4.3 Assumptions

N/A

#### 14.4.4 Access Routine Semantics

initiateScan():

- transition:
  - Receives the input 'img' from the user to initiate the scanning process.
- output: N/A
- exception:

- Throws 'InvalidImageFormatException' if a non-JPEG/JPG image is used as input.

#### **14.4.5 Local Functions**

N/A

## 15 MIS of View Results Module

### 15.1 Module

ViewResults

### 15.2 Uses

N/A

### 15.3 Syntax

#### 15.3.1 Exported Constants

N/A

#### 15.3.2 Exported Access Programs

Name	In	Out	Exceptions
displayReport	report: Radiology Report	-	-

### 15.4 Semantics

#### 15.4.1 State Variables

N/A

#### 15.4.2 Environment Variables

N/A

#### 15.4.3 Assumptions

N/A

#### 15.4.4 Access Routine Semantics

displayReport():

- transition:
  - Displays the generated radiology report on the GUI.
- output: N/A
- exception: N/A

### 15.4.5 Local Functions

N/A

## 16 MIS of AI Model Module

### 16.1 Module

AIModel

### 16.2 Uses

- ChestXRayRead
- ResultGen

### 16.3 Syntax

#### 16.3.1 Exported Constants

N/A

#### 16.3.2 Exported Access Programs

Name	In	Out	Exceptions
processImg	img: JPEG image	procImg: Processed Image	InvalidImageFormatException
generateResults	procImg: Processed Image	classification: Disease Classification	-

### 16.4 Semantics

#### 16.4.1 State Variables

N/A

#### 16.4.2 Environment Variables

N/A

#### 16.4.3 Assumptions

N/A



#### 16.4.4 Access Routine Semantics

processImg():

- transition:
  - Uses the trained model to process the given 'img'.
- output:
  - 'procImg' contains the processed information and findings from the chest X-ray analysis.
- exception:
  - Throws 'InvalidImageFormatException' if the provided 'image' is not a JPEG or JPG image.

generateResults():

- transition:
  - Uses the trained model to interpret the processed image and generate a disease classification.
- output:
  - 'classification' contains the generated disease classification for each disease.
- exception: N/A

#### 16.4.5 Local Functions

N/A

## 17 MIS of NLP Model Module

### 17.1 Module

NLPModel

### 17.2 Uses

- RepCompGen

### 17.3 Syntax

#### 17.3.1 Exported Constants

N/A

#### 17.3.2 Exported Access Programs

Name	In	Out	Exceptions
generateReport	report: Radiology Report	nlp: NLP Report	InvalidReportFormatException

### 17.4 Semantics

#### 17.4.1 State Variables

N/A

#### 17.4.2 Environment Variables

N/A

#### 17.4.3 Assumptions

N/A

#### 17.4.4 Access Routine Semantics

generateReport():

- transition:
  - Uses the RepCompGen module to generate a radiology report.
- output:

- 'nlp' contains the generated NLP report.
- exception:
  - Throws 'InvalidReportFormatException' if the provided 'report' is in an invalid format.

#### **17.4.5 Local Functions**

N/A

## 18 MIS of Backend Module

### 18.1 Module

Backend

### 18.2 Uses

- UserAuthMgmt
- MedInstInter
- DatabaseOps

### 18.3 Syntax

#### 18.3.1 Exported Constants

N/A

#### 18.3.2 Exported Access Programs

Name	In	Out	Exceptions
connectDatabase	credentials: str	connectionStatus: bool	InvalidCredentialsException
disconnectDatabase		success: bool	-

### 18.4 Semantics

#### 18.4.1 State Variables

#### 18.4.2 Environment Variables

#### 18.4.3 Assumptions

#### 18.4.4 Access Routine Semantics

connectDatabase():

- transition: N/A
- output:
  - 'connectionStatus' is set to True if the connection is successful, False otherwise.
- exception:

- Throws 'InvalidCredentialsException' if the provided credentials are invalid.

disconnectDatabase():

- transition: N/A
- output:
  - 'success' is set to True if the disconnection is successful, False otherwise.
- exception: N/A

#### **18.4.5 Local Functions**

N/A

## 19 MIS of App Controller Module

### 19.1 Module

AppController

### 19.2 Uses

- AIModel
- NLPModel
- AppGUI
- Backend

### 19.3 Syntax

#### 19.3.1 Exported Constants

N/A

#### 19.3.2 Exported Access Programs

Name	In	Out	Exceptions
accessBackend	-	-	-
accessGUI	-	-	-
accessAI	-	-	-
accessNLP	-	-	-

### 19.4 Semantics

#### 19.4.1 State Variables

N/A

#### 19.4.2 Environment Variables

N/A

#### 19.4.3 Assumptions

N/A

#### 19.4.4 Access Routine Semantics

accessBackend():

- transition:
  - Controller accesses the backend server.
- output: N/A
- exception: N/A

accessGUI():

- transition:
  - Controller accesses the application GUI.
- output: N/A
- exception: N/A

accessAI():

- transition:
  - Controller accesses the AI Model.
- output: N/A
- exception: N/A

accessNLP():

- transition:
  - Controller acceses the NLP Model.
- output: N/A
- exception: N/A

#### 19.4.5 Local Functions

N/A

## References

- Carlo Ghezzi, Mehdi Jazayeri, and Dino Mandrioli. *Fundamentals of Software Engineering*. Prentice Hall, Upper Saddle River, NJ, USA, 2nd edition, 2003.
- Daniel M. Hoffman and Paul A. Strooper. *Software Design, Automated Testing, and Maintenance: A Practical Approach*. International Thomson Computer Press, New York, NY, USA, 1995. URL <http://citeseer.ist.psu.edu/428727.html>.



## 20 Appendix